



Hardware Triggered Scanning: Scanning Components

Philip Taylor, Emma Arandjelovic
Observatory Sciences Limited



Overview


- To perform scans, we need to control complex devices such as motion controllers and detectors
- This topic will cover how these devices are configured by Malcolm
- We will use a simulated motion controller and detector to see how this works in practice
- We will also look at how to generate scan trajectories



Block Structure


- Blocks like we have created so far make up the lowest level of blocks, known as the *Hardware Layer*
- The *Device Layer* is a higher level of blocks that synchronise a number of child blocks to create an interface to a device e.g. a detector or motor controller

Example:



Device Blocks

- Device blocks instantiate one or more child blocks
- A controller aimed at managing child blocks is called a *ManagerController*
- This provides:
 - A *layout* of child blocks
 - A mechanism to enable / disable child blocks
 - A mechanism to load / save designs (saved configurations)
- State machine:






Motion Demo

- Two *Counter* blocks in the hardware layer simulate motor axes (X and Y)
- A device block called *Motion* simulates a motor controller with *XMove* and *YMove* methods
- The *Counter* blocks are constructed with *BasicControllers*
- The *Motion* block is constructed with a *ManagerController*

See the Motion tutorial:

<https://pymalcolm.readthedocs.io/en/latest/tutorials/motion.html>





Motion Block Definition

demo/blocks/motion_block.yaml (extract):


```
- builtin.controllers.ManagerController:
  mri: $(mri)
  config_dir: $(config_dir)
  description: $(docstring)

- demo.blocks.counter_block:
  mri: $(mri):COUNTERX
- demo.blocks.counter_block:
  mri: $(mri):COUNTRY

- demo.parts.CounterMovePart:
  name: x
  mri: $(mri):COUNTERX
- demo.parts.CounterMovePart:
  name: y
  mri: $(mri): COUNTRY
```

- Use a *ManagerController* to construct the block
 - This takes an additional parameter to specify a directory for saving configurations
- Instantiate two *Counter* blocks
- Instantiate two *MotorMove* parts, one per block


Motion Example Structure



Notice the *design* attribute and *save* method are contributed by the *ManagerController*

Designs

Recall: Designs allow you to save the layout of a device block and the attributes of its child blocks



The screenshot displays the PANDABOX software interface. The main workspace shows a design layout for 'BL49P-ML-SCAN-01' with three device blocks: DET (GigE Vision mri), PANDABOX (PandA), and BRICK (Brick with X and A Sample stage motors). The right-hand panel shows the properties for 'BL49P-ML-PANDA-01'. A red circle highlights the 'Design' dropdown menu, which lists various design templates: 2seq_pcomp, Mapping Test #1, edge_triggered_pcap, level_triggered_pcap, pmac_master, template_double_seq_pcomp, and template_live_dead_framed_pcap. The 'Design' dropdown is currently set to 'PandA'.



Designs continued

- Which attributes are saved?
- Recall the attribute definition in our *Counter* part:

```
self.counter = NumberMeta(  
    "float64", "The current value of the counter",  
    tags=[config_tag(), Widget.TEXTINPUT.tag()])
```

- Tags contain important metadata:
 - Widget tags - choose which widget to display it with
 - Group tags - create an expander around related attributes
 - Port tags - create the visual connections between blocks
 - Config tags - turn on saving/loading

See: <https://pymalcolm.readthedocs.io/en/latest/reference/tags.html>



How to Create a Device Part

1. Subclass *builtin.parts.ChildPart*
 - This provides access to the child blocks
2. Create annotypes for any variables we need
3. Create an *__init__()* method:
 - Call *super().__init__()*
4. Create a *setup()* method:
 - Call *super().setup()*
 - Register any additional methods and attributes
5. Implement the methods to perform the business logic
 - Use the passed in *context* to access the child block



Motor Move Part Definition

demo/parts/countermovepart.py (extract):

```
with Anno("The demand value to move our
counter motor to"):
    ADemand = float

@builtin.util.no_save("counter")
class CounterMovePart(builtin.parts.ChildPart):
    def __init__(self, name, mri):
        # type: (APartName, builtin.parts.AMri) -> None
        super(CounterMovePart, self).__init__(name,
            mri, initial_visibility=True, stateful=False)
    def setup(self, registrar):
        # type: (PartRegistrar) -> None
        super(CounterMovePart, self).setup(registrar)
        registrar.add_method_model(self.move,
            self.name + "Move", needs_context=True)
```

- Create an Annotype for the demand position
- Define attributes in the child block not to be saved
- Call `super().__init__()`
- Register the method model in **setup**, ensuring the name is unique for each part
- Tell Malcolm the method needs a *Context* object for accessing the child block



Part Definition Continued

demo/parts/countermovepart.py (continued):

```
@add_call_types
def move(self, context, demand):
    # type: (builtin.hooks.AContext, ADemand) -> None
    child = context.block_view(self.mri)
    child.counter.put_value(demand)
```

1. Use the passed in context as well as the MRI to create a Block view of the *Counter* child Block
2. Set the Block's *counter* attribute to the new demand




Running the Motion Demo



`./malcolm/imalcolm.py`

`malcolm/modules/demo/DEMO-MOTION.yaml`


1. Open <http://localhost:8008/>
2. Select the MOTION:COUNTERX block and 'tear-off'
3. Repeat for MOTION:COUNTRY
4. Select the MOTION root block and request new values for X and Y
5. Watch the counter blocks changing value
6. Practise saving and loading designs with different *delta* values





Detector Demo

- A dummy detector that takes a scan specification and writes data to an HDF5 file



- This will introduce two new concepts:
 1. *RunnableController* – provides methods and attributes relating to the different steps of a scan
 2. *ScanPointGenerator* – used to generate n-dimensional scan paths

See the Detector tutorial:

<https://pymalcolm.readthedocs.io/en/latest/tutorials/detector.html>



Runnable Device Blocks

- A 'Runnable' device block understands what a scan is
- It has a state machine that implements the scan process
- This is controlled by the user via two important methods:
 - `configure()`
Configure all child blocks according to the supplied parameters
 - `run()`
Start all child blocks running, provide status monitoring and periodic actions
- These are contributed to the block by a *RunnableController*



Runnable Controller Details

- Inherits from *ManagerController*
- Adds to the *ManagerController* state machine to include the different steps of a scan
- Additional attributes:
 - *completedSteps*
 - *configuredSteps*
 - *totalSteps*
- Controllers provide *hooks*, which Parts can register with to run the device specific logic on state transitions




Runnable Controller States

Hooks are run on state transitions

See:

<https://pymalcolm.readthedocs.io/en/latest/reference/state-sets.html#state-sets>






Hooks Overview

- Functions registered to hooks are called automatically by the controller
- All functions hooked to the same hook are called concurrently
- Examples:
 - **ConfigureHook** (part, context, completed_steps, steps_to_do, part_info, generator, axesToMove, **kwargs):
 - Called at configure() to setup a child block
 - **RunHook** (part, context, **kwargs):
 - Called at run() to start all children running
 - **ResetHook**(part)
 - Called at reset() to return to Ready state



Runnable Device Sequence Diagrams




Note: There are more hooks than are shown here, see <https://pymalcolm.readthedocs.io/en/latest/reference/statesets.html> for a complete list of hooks.



Detector Demo Structure

- The *DETECTOR* block is constructed with a *RunnableController*
- A *FileWritePart* writes the HDF5 file
 - *Configure hook takes a scan point generator and initialises the scan*
 - *Run hook iterates through the generator points and writes the data*
- A *DatasetTablePart* collects info (location, size etc.) about the data and exposes it in the *datasets* attribute





Detector Block Definition


demo/blocks/detector_block.yaml (extract):

```
- builtin.controllers.RunnableController:  
  mri: $(mri)  
  config_dir: $(config_dir)  
  description: $(docstring)  
  
- builtin.parts.LabelPart:  
  value: $(label)  
  
- scanning.parts.DatasetTablePart:  
  name: DSET  
  
- demo.parts.FileWritePart:  
  name: FW  
  width: 160  
  height: 120
```

- Use a *RunnableController* to construct the block
- A *LabelPart* displays a title on the GUI
- A *DatasetTablePart* reports the datasets written
- A *FileWritePart* writes some dummy data to an HDF5 file



Detector Demo Sequence Diagram



This demo detector saves some generated data to an HDF5 file just as an example.



How to Create a Runnable Device Part

1. Subclass *builtin.parts.Part* or *builtin.parts.ChildPart*
2. Create an `__init__` method:
 - Call *super().__init__()*
3. Implement the *on_configure(*params)*, *on_run()* and *reset()* methods to perform the business logic
 - The *on_configure* method should prepare the scan
 - The *on_run* method should generate the data and report progress to the controller so it can update the block's *currentStep* attribute
 - The *reset* method should clean up after a scan, even if it failed.
4. Create a *setup* method:
 - Register any additional methods and attributes
 - Register any hooks with the controller
 - Inform the controller what extra parameters *configure* requires




File Write Part: init

demo/parts/filewritepart.py (extract):

```
class FileWritePart(Part):  
    def __init__(self, name, width, height):  
        # type: (APartName, AWidth, AHeight) -> None  
        super(FileWritePart, self).__init__(name)  
        # Store input arguments  
        self._width = width  
        self._height = height  
        # The detector image we will modify for each image (0..255 range)  
        self._blob = make_gaussian_blob(width, height) * 255  
        # The HDF5 file we will write  
        self._hdf = None          # type: h5py.File  
        # Configure args and progress info  
        self._generator = None    # type: scanning.hooks.AGenerator  
        self._completed_steps = 0  
        self._steps_to_do = 0  
        # How much to offset unique id value from generator point (for rewind function)  
        self._uid_offset = 0
```




Detector Demo Sequence Diagram



This demo detector saves some generated data to an HDF5 file just as an example.




File Write Part: on_configure

demo/parts/filewritepart.py (extract):

```
def on_configure(self,
    completed_steps,      # type: scanning.hooks.ACompletedSteps
    steps_to_do,          # type: scanning.hooks.AStepsToDo
    generator,            # type: scanning.hooks.AGenerator
    fileDir,              # type: scanning.hooks.AFileDir
    formatName="det",     # type: scanning.hooks.AFormatName
    fileTemplate="%s.h5", # type: scanning.hooks.AFileTemplate
):
    # type: (...) -> scanning.hooks.UInfos
    # Store the parameters and create an empty HDF file
    self._generator = generator
    filename = fileTemplate % formatName
    filepath = os.path.join(fileDir, filename)
    self._hdf = self._create_hdf(filepath, generator)
    # Create the Info objects describing the two datasets.....
    infos = [scanning.infos.DatasetProducedInfo(name=...) , ...]
    return infos
```



Detector Demo Sequence Diagram



This demo detector saves some generated data to an HDF5 file just as an example.




File Write Part: on_run

demo/parts/filewritepart.py (extract):

```
def on_run(self, context):
    # type: (scanning.hooks.AContext) -> None
    """On `RunHook` record where to next take data"""
    end_of_exposure = time.time() + self._exposure # Start time so everything is relative
    last_flush = end_of_exposure
    for i in range(self._completed_steps, self._completed_steps + self._steps_to_do):
        point = self._generator.get_point(i)
        wait_time = end_of_exposure - time.time()
        context.sleep(wait_time) # Must use the context parameter to make it
                                # an interruptible sleep so the scan can be aborted
        self.log.debug("Writing data for point %s", i)
        self._write_data(point, i)
        if time.time() - last_flush > FLUSH_PERIOD:
            last_flush = time.time()
            self._flush_datasets()
            end_of_exposure += point.duration
        self.registrar.report(scanning.infos.RunProgressInfo(i + 1))
    self._flush_datasets()
```



Detector Demo Sequence Diagram



This demo detector saves some generated data to an HDF5 file just as an example.



File Write Part: reset

demo/parts/filewritepart.py (extract):

```
def on_reset(self):  
    # type: () -> None  
    if self._hdf:  
        self._hdf.close()  
        self._hdf = None
```



File Write Part: setup

demo/parts/filewritepart.py (continued):

```
def setup(self, registrar):
    # type: (PartRegistrar) -> None
    super(FileWritePart, self).setup(registrar)
    # Hooks
    registrar.hook(scanning.hooks.ConfigureHook, self.on_configure)
    registrar.hook((scanning.hooks.PostRunArmedHook,
                    scanning.hooks.SeekHook), self.on_seek)
    registrar.hook(scanning.hooks.RunHook, self.on_run)
    registrar.hook((scanning.hooks.AbortHook,
                    builtin.hooks.ResetHook), self.on_reset)
    # Tell the controller to expose some extra configure parameters
    registrar.report(scanning.hooks.ConfigureHook.create_info(self.on_configure))
```



Additional Hooks

More examples:

- **SeekHook**
 - Used for pause/rewind functionality, to adjust the *completed_steps*
- **PostRunArmedHook**
 - Called at the end of run() when there are more steps to be run (e.g. when the scan is paused)
- **AbortHook**
 - Called at abort() to stop the current scan
- **ResetHook**
 - Called at reset() to reset all Parts to a known state



Info Objects

- An *Info* is an object created by a Part and passed to its parent Controller in one of two ways:
 - Returned from a hooked function
 - Reported via *PartRegistrar.report(info)*
- We report progress via *PartRegistrar.report(info)* in the method registered to the RunHook:

```
def on_run(self, context):  
    ...  
    for i in range(...):  
        ...  
        self.registrar.report(  
            scanning.infos.RunProgressInfo(i+1))
```



Info Objects Continued

- We return a description of the dataset produced by the *Part* from the method registered to the *ConfigureHook*:

```
def on_configure(self, ...):  
    # type: (...) -> scanning.hooks.UInfos  
    ...  
    infos = [scanning.infos.DatasetProducedInfo(name=...), ...]  
    return infos
```

- *DatasetProducedInfo* parameters:
 - **name** – Dataset name
 - **filename** – Filename relative to the fileDir we were given
 - **type** – What NeXuS dataset type it produces
 - **rank** – The rank of the dataset including generator dims
 - **path** – The path of the dataset within the file
 - **uniqueid** – The path of the UniqueID dataset within the file




Info Objects Continued

- We report additional arguments for hook functions in the *setup* method:
def setup(self, registrar):
 ...
 registrar.report(
 scanning.hooks.**ConfigureHook.create_info**(self.configure))
- *create_info* scans the *on_configure* method for the arguments in addition to *completed_steps* and *steps_to_do*:
 - *generator*
 - *fileDir*
 - *formatName*
 - *fileTemplate*



Detector Demo

- A dummy detector that takes a scan specification and writes data to an HDF5 file



- This will introduce two new concepts:
 1. *RunnableController* – provides methods and attributes relating to the different steps of a scan
 2. *ScanPointGenerator* – used to generate n-dimensional scan paths

See the Detector tutorial:

<https://pymalcolm.readthedocs.io/en/latest/tutorials/detector.html>



Scan Point Generator

- To run a continuous scan, we need many parameters:
 - Demand positions for frame midpoint
 - Demand positions for frame start and end
 - Index at which to store the frame
 - Duration of the frame
- A Python module called *ScanPointGenerator* generates this data for a variety of scan types (e.g. line, spiral, ...)
- A generator is passed to the *configure* method to setup the scan
- See <https://scanpointgenerator.readthedocs.io>



Scan Point Generator Example

class scanpointgenerator.LineGenerator

(axes, units, start, stop, size, alternate=False)

- Generate a line of equally spaced N-dimensional points

Parameters:

- **axes (str/list(str))** – The scannable axes E.g. “x”
- **units (str/list(str))** – The scannable units. E.g. “mm”
- **start (float/list(float))** – The first position to be generated. e.g. 1.0
- **stop (float/list(float))** – The final position to be generated. e.g. 5.0
- **size (int)** – The number of points to generate. E.g. 5
- **alternate (bool)** – Specifier to reverse direction



Preparing the Example



`./malcolm/imalcolm.py`
`malcolm/modules/demo/DEMO-DETECTOR.yaml`

```
>>> from scanpointgenerator import LineGenerator, CompoundGenerator
>>> from scanpointgenerator.plotgenerator import plot_generator
>>> yline = LineGenerator("y", "mm", -1, 0, 6)
>>> xline = LineGenerator("x", "mm", 4, 5, 5, alternate=True)
>>> generator = CompoundGenerator([yline, xline], duration=0.5)
>>> plot_generator(generator)
>>> from annotypes import json_encode
>>> json_encode(generator)
```


- Copy the JSON output to the clipboard



Running the Example



- Open <http://localhost:8008/gui/DETECTOR>
- Expand the *Configure* method and *Edit* the generator field to paste in the JSON code
- Set the *fileDir* to */tmp*
- Click *Configure*
- Click *Run*
- Watch the *Completed Steps counter* increase






Checking the Results



```
>>> import h5py
>>> with h5py.File("/tmp/det.h5") as f:
.....:     im = f[("/entry/sum")][:]
.....:
>>> im.shape
(6, 5, 1, 1)
>>> from pylab import imshow, show
>>> imshow(im[:, :, 0, 0])
<matplotlib.image.AxesImage at 0x7f11b090>
>>> show()
```





Recap: yaml chains

```
DEMO-DETECTOR.yaml x
1 # Create some Blocks
2 - demo.blocks.detector_block:
3   mri: DETECTOR
4   config_dir: /tmp
5
6 # Add a webserver
7 - web.blocks.web_server_block:
8   mri: WEB
```

```
DEMO-DETECTOR.yaml x *detector_block.yaml x
1 - builtin.parameters.string:
2   name: mri
3   description: MRI for created block
4
5 - builtin.parameters.string:
6   name: config_dir
7   description: Where to store saved configs
8
9 - builtin.parameters.int32:
10  name: width
11  description: Width of the produced image
12  default: 160
13
14 - builtin.parameters.int32:
15  name: height
16  description: Height of the produced image
17  default: 120
18
19 - demo.parts.FileWritePart:
20  name: FW
21  width: $(width)
22  height: $(height)
23
```

```
DEMO-DETECTOR.yaml x detector_block.yaml x filewritepart.py x
28
29 class FileWritePart(Part):
30     """Minimal interface demonstrating a file writing
31     detector part"""
32     def __init__(self, name, width, height):
33         # type: (APartName, AWidth, AHeight) -> None
34         super(FileWritePart, self).__init__(name)
35         # Store input arguments
36         self._width = width
37         self._height = height
38         # The detector image we will modify for each image
39         # (0..255 range)
40         self._blob = make_gaussian_blob(width, height) * 255
41         # The hdf file we will write
42         self._hdf = None # type: h5py.File
43         # Configure args and progress info
44         self._exposure = None
45         self._generator = None # type:
46         scanning.hooks.AGenerator
47         self._completed_steps = 0
48         self._steps_to_do = 0
49         # How much to offset uid value from generator point
50         self._uid_offset = 0
```




Recap: yaml chains

```
DEMO-DETECTOR.yaml x
1 # Create some Blocks
2 - demo.blocks.detector_block:
3   mri: DETECTOR
4   config_dir: /tmp
5
6 # Add a webserver
7 - web.blocks.web_server_block:
8   mri: WEB
```

```
DEMO-DETECTOR.yaml x *detector_block.yaml x
1 - builtin.parameters.string:
2   name: mri
3   description: MRI for created block
4
5 - builtin.parameters.string:
6   name: config_dir
7   description: Where to store saved configs
8
9 - builtin.parameters.int32:
10  name: width
11  description: Width of the produced image
12  default: 160
13
14 - builtin.parameters.int32:
15  name: height
16  description: Height of the produced image
17  default: 120
18
19 - demo.parts.FileWritePart:
20  name: FW
21  width: $(width)
22  height: $(height)
```

```
DEMO-DETECTOR.yaml x detector_block.yaml x filewritepart.py x
28
29 class FileWritePart(Part):
30     """Minimal interface demonstrating a file writing
31     detector part"""
32     def __init__(self, name, width, height):
33         # type: (APartName, AWidth, AHeight) -> None
34         super(FileWritePart, self).__init__(name)
35         # Store input arguments
36         self._width = width
37         self._height = height
38         # The detector image we will modify for each image
39         # (0..255 range)
40         self._blob = make_gaussian_blob(width, height) * 255
41         # The hdf file we will write
42         self._hdf = None # type: h5py.File
43         # Configure args and progress info
44         self._exposure = None
45         self._generator = None # type:
46         # Scanning hooks: AGenerator
47
```





Practical Exercises



4. Modify the Motion demo to add high and low soft limits (values to be provided in DEMO-MOTION.yaml). If a demand is requested out of this range, it should raise an error. Provide default values of 0, meaning soft limits are disabled.
Hint 1: Start by modifying the DEMO-MOTION.yaml file and work your way down.
Hint 2: Read the docs! See builtin parameter types from second page:
<https://pymalcolm.readthedocs.io/en/latest/tutorials/motion.html>
https://pymalcolm.readthedocs.io/en/latest/build/builtin/parameters_a_pi.html
5. Modify the generator in the Detector demo to create a larger and faster scan.