



Hardware Triggered Scanning: Hardware


Philip Taylor, Emma Arandjelovic
Observatory Sciences Limited

Overview

Motion trajectory
control

Flexible triggering, and
fast position capture


Data capture





Typical Triggering Setup

Trigger at start of each frame




Detector exposing

Detector readout



Triggering Setup: Warning

Trigger received when
detector still reading out...






Triggering Setup: Warning

Trigger received when
detector still reading out...

... leads to missing frames!!






Trigger Signals

- The PMAC has 3 GPIO **output** signals, generated during a scan
- They are fed into the PandA box, which in turn triggers detector acquisition
- There are a wide range of detectors in use at Diamond which may require different trigger signals – the PandA can be configured to output an appropriate signal
- The GPIO signals are:
 1. *Live*: The start of a frame
 2. *Dead*: The start of a period when no frames are acquired
 3. *Centre*: The middle of the current frame (not normally used)

Snake Scan





Motion Control Hardware

- Motion control performed by Delta Tau PMAC
- Note: Must be CLIPPER or Geobrick LV IMS II
(old style Geobricks don't have enough memory)
- Co-ordinate system needed to run motion program





Trajectory Control Overview

- Precise control of the trajectory is needed
- For each step, the user specifies the end position, the velocity at the end point and the time period for that step
- This is called **PVT** mode:
 - **P**osition
 - **V**elocity
 - **T**ime



Trajectory Control Sequence


1. Malcolm generates trajectory positions and delta times using the Python module *scanpointgenerator*
2. EPICS writes them to the PMAC via the *pmac* driver
3. The PMAC calculates velocities and performs the scan using PVT mode
4. PMAC User Subroutines, which output GPIO signals, are called as each scan point is reached. 8 routines are available, corresponding to all combinations of the 3 signals: Live + Dead + Centre




PandABox Overview

- *Position* and *Acquisition* Box
- Flexible, modular platform
- Supports a variety of scanning applications
- Based on a PicoZed module which integrates:
 1. ARM processor running embedded Linux
 2. FPGA firmware providing configurable function blocks

PandA Hardware



Functional Blocks



- Each one performs a specific function, e.g. position capture
- Configured via a number of registers
- Wired up and configured at runtime via a TCP server
- Configurations saved to disk



Available Function Blocks

| Function Block | Description |
|----------------|--|
| INENC, OUTENC | Input and output encoder signals |
| TTLIN, TTLOUT | TTL Input and output signals |
| COUNTER | Up/down pulse counter |
| CALC | Sums up to 4 position inputs |
| LUT | Performs logical function on 5 inputs |
| PCAP | Position capture – options for min, max, sum |
| PULSE | One-shot pulse delay and stretch |
| SRGATE | Set Reset Gate – set configurable high or low output |
| | |

For the full list and more information, see:

<https://pandablocks-fpga.readthedocs.io/en/latest/blocks.html>



Configuring PandA

- PandA runs an embedded web server
- Point a browser directly to this for low level configuration such as updating packages
- For normal use, use Malcolm web gui:
 - Test rig: <http://localhost:8008/gui>
 - Beamlines: <http://ixx-control:8008/gui>

