



Hardware Triggered Scanning: Troubleshooting Solutions

Philip Taylor, Emma Arandjelovic Observatory Sciences Limited



Scenario 1

Symptoms

Scan fails immediately with message:

Error from Malcolm Device Connection: Failed to connect to device 'BL47P-ML-SCAN-01' (ERROR: channel not connected)

Troubleshooting techniques

- Check the message from GDA this indicates a problem with Malcolm
- Check the Malcolm web GUI is responding – http://localhost:8008/gui

Problem: Malcolm is not running

Solution

Restart the Malcolm process:

ioc-connect BL4xP-ML-MALC-01 Ctrl-X



Scenario 2

Symptoms

Scan running extremely slowly, possibly eventually failing with a TimeoutError

Troubleshooting techniques

- Check the motors are moving at the correct speed
- Check the frames are being written to disk (detector EDM screen)


Problem: File writing couldn't keep up with the frame rate

Solution

- Check the HDF chunking settings are valid
- If in doubt, set the column and row chunk sizes to 0 (defaults to frame size)



Scenario 2 Screenshot





Scenario 3

Symptoms

Scan fails immediately with a message something like:

org.eclipse.scanning.api.malcolm.MalcolmDeviceException: Error from Malcolm Device Connection: NotWriteableError: Field ['BL49P-ML-SCAN-01', 'reset'] is not writeable, maybe because Block state = Aborting, cause:

Troubleshooting techniques

- Check the message from GDA this indicates problem is downstream from Malcolm
- Check Malcolm GUI shows PVs disconnected

Problem: EPICS IOC is not running

Solution


Reboot the EPICS IOC:


ioc-connect BL4xP-EA-IOC-01 Ctrl-X



Scenario 3 Screenshots

Check the status of some EPICS devices in the Malcolm web gui







Scenario 4 (I)

Symptoms

Scan fails with error message:

org.eclipse.scanning.api.malcolm.MalcolmDeviceException: Error from Malcolm Device Connection: AssertionError: Velocity 499.750124938 invalid for 'stagea' with max_velocity 10.0, cause:

Troubleshooting techniques

 Check the motor record for stagea is setup correctly

Problem: VMAX is too small.

Solution

Set VMAX for stagea back to its original value (look in /tmp/ex4_values.txt)

From Training Rig edm screen: Motors -> Under A (right hand side) More -> Motion



Scenario 4 (II)

Symptoms

Scan is very slow and possibly times out

Troubleshooting techniques

- Notice that the stagea motor rotates and then everything stops -> something wrong with the stagex
- Check the VMAX motor record for stagex

Problem: VMAX too small but this time for stagex

Solution

Same as before. Correct value is in /tmp/ex4_values.txt.

Notice that the same underlying problem has different symptoms for each axis. Can you explain why?

Please run

> sh /tmp/ex4_restore.sh

1/6/2020