



Hardware Triggered Scanning: Course 2

Philip Taylor, Emma Arandjelovic
Observatory Sciences Limited



Course Aims

1. Understand the key concepts behind Malcolm
2. Get experience in setting up Malcolm from scratch
3. Understand the different strategies for hardware triggered scanning on beamlines
4. Learn how to configure Malcolm to run on beamlines



Course Content

1. Refresher and basic concepts
 - Blocks, parts, and controllers
2. Scanning components
 - Manager and Runnable controllers
 - Scan Point Generator
 - Designs
3. Scanning control
 - Scan layer
 - AreaDetector
4. From the classroom to the beamline
 - Scanning strategies
 - EPICS and motion control requirements
 - Configuring Malcolm to run on beamlines



Refresher: Software Stack




Data Analysis WorkbeNch
- Analysis and visualization



Generic Data Acquisition
- Experiment setup and supervision




Malcolm
- Scan configuration



Experimental Physics &
Industrial Control System
- Low level control of hardware




Refresher: What is Malcolm?

- Generic and extensible framework for scanning
- Middle layer between GDA and control system
- Implemented in  python™
- Creates a s/w map of the h/w layer
- <https://pymalcolm.readthedocs.io>
- Web GUI called MalcolmJS
- <https://malcolmjs.readthedocs.io>



Refresher: Blocks



- A **block** is a user-centred view
Examples:
 - Motion controller
 - PandA
 - Detector
 - 'Hello World' program
- It is an *interface* to an object's:
 - **Methods** (actions)
 - **Attributes** (data)



Blocks Continued...


- A process hosts multiple blocks
- Blocks are given a unique **MRI**
(Malcolm Resource Identifier)
 - This is used by clients to address the block
- A block is just the interface
 - Contains no code!



Parts and Controllers


- Methods and attributes are contributed in Python **Parts**
- A **Controller** provides a co-ordinating framework:

- It creates a **Block View** that we interact with
- It populates this with its own methods and attributes and also those from the contributed Parts



Stateful Controller

- Used by blocks in the Hardware Layer
- Implements a state machine:



Unlabelled transitions take place in response to internal actions

Labelled transitions are triggered externally.



Getting Started



- Stop BL4xP-ML-MALC-01 (Ctrl-T then Ctrl-X)
- Get the code:

git clone <https://github.com/dls-controls/pymalcolm>

- Run the Hello World example:

`./malcolm/imalcolm.py malcolm/modules/demo/DEMO-HELLO.yaml`

NB: The examples are based on tutorials which you can also try in your own time:

<https://pymalcolm.readthedocs.io>



Demo 1: Hello World



imalcolm launches an *IPython* interactive console

In the console, try:

```
>>> hello = self.block_view("HELLO")
```

```
>>> hello.greet("me")
```

```
Manufacturing greeting...
```

```
'Hello me'
```

```
>>> result = hello.greet("Emma")
```

```
Manufacturing greeting...
```


```
>>> print(result)
```

```
Hello Emma
```

- *self* is a context we can use to get hold of block views
- “*HELLO*” is the MRI of our block
- The *hello* object is now a view of our HELLO block that we can use to call its methods
- The return value is assigned to the variable *result*



Module Structure



- Each item identified using the ‘.’ separator, e.g.
demo.blocks.hello_block and *demo.parts.HelloPart*



Looking Under the Hood


DEMO-HELLO.yaml

Create some Blocks

- demo.blocks.hello_block:
 mri: HELLO
- demo.blocks.hello_block:
 mri: HELLO2
- demo.blocks.counter_block:
 mri: COUNTER

Add a webserver

- web.blocks.web_server_block:
 mri: WEB



Instantiates 3 blocks. Each one defines a unique **Malcolm Resource Identifier (MRI)**

Web server block starts an HTTP server on port 8008



Hello Block Definition


demo/blocks/hello_block.yaml:

```
- builtin.parameters.string:
  name: mri
  description: Malcolm resource id

- builtin.defines.docstring:
  value: Block with a greet() Method

- builtin.controllers.BasicController:
  mri: $(mri)
  description: $(docstring)

- demo.parts.HelloPart:
  name: hello
```



Defines parameters to be defined when instantiating the block. Value obtained using `$(mri)`

Defines the `$(docstring)` variable which describes the block's function


Defines the controller which will create the block for us. The *BasicController* is a simple container for parts

Defines the part which contributes the functionality of the block

- *name* is the (unique) name of the part within the controller



Hello World Structure



- The BasicController contributes a 'health' attribute
- The business logic methods are contributed by our HelloPart



Brief Aside: Annotypes

- Library for annotating Python types with metadata
- Allow clients to discover parameter type information at runtime

Info	
Parameter Type	Input
typeid	malcolm:core/Number Meta:1.0
Description	Time to wait before returning
Required?	false
Default Value	0
DISCARD PARAM	



Brief Aside: Annotypes

In the Python part:

```
from annotypes import Anno, add_call_types  
with Anno("my variable"):
```

```
    AMyVar = str
```

Variable name

Variable type

Variable description

add_call_types is a Python decorator which uses our Annotype definitions to provide introspection information



Brief Aside: Annotypes

with Anno("Description of parameter 1 (a string)":

AParam1 = str

with Anno("Description of parameter 2 (an int)":

AParam2 = int

with Anno("Description of the return value (a float)":

AReturn = float

@add_call_types

def **myMethod**(self, param1, param2):

type: (AParam1, AParam2) -> AReturn

See: <https://github.com/dls-controls/annotypes>




Hello Part Definition

demo/parts/hellopart.py:

```
with Anno("The name of the person to greet"):
    AName = str
with Anno("Time to wait before returning"):
    ASleep = float
with Anno("The manufactured greeting"):
    AGreeting = str

class HelloPart(Part):
    def setup(self, registrar):
        # type: (PartRegistrar) -> None
        super(HelloPart, self).setup(registrar)
        registrar.add_method_model(self.greet)
        registrar.add_method_model(self.error)
```



- Type definition information
- Class must extend *Part* and provide a *setup* method
- The *PartRegistrar* is a utility object used to register methods and attributes with the block



Hello Part Definition

demo/parts/hellopart.py (continued...):

```
@add_call_types
def greet(self, name, sleep=0):
    # type: (AName, ASleep) -> AGreeting
    """Optionally sleep <sleep> seconds, then
       return a greeting to <name>"""
    print("Manufacturing greeting...")
    sleep_for(sleep)
    greeting = "Hello %s" % name
    return greeting

def error(self):
    """Raise an error"""
    raise RuntimeError("You called error()")
```

- *greet* method contains the 'business logic' of the part
- *error* method can be used to raise an error if something goes wrong
- The three Annotypes we defined earlier are used to provide introspection information



Connecting a Client



Start up a second process to connect to the first:


```
./malcolm/imalcolm.py -c ws://localhost:8008
```

```
>>> self.make_proxy("localhost:8008", "HELLO")  
>>> self.block_view("HELLO").greet("me")  
u'Hello me'
```

Check the output from the first process to see this is the one doing the actual “work”



Connecting a Client






Demo 2: Adding Attributes

- We have seen how to add a method to a block
 - registrar.add_method_model(<method>)
- To add attributes we need a little more information:
 - registrar.add_attribute_model(
 "<attribute_name>",
 <attribute_model_instance>, —————> Contains type info for validation
 <method> —————> Method called when
 someone tries to set the value

Counter Demo Structure

- a writeable attribute `counter` which keeps the current counter value
- a writeable attribute `delta` which stores the amount to increment by
- a method `zero()` which will set `counter = 0`
- a method `increment()` which will set `counter = counter + delta`




➤ Notice the additional attribute 'health' contributed by the BasicController



Counter Part Definition

demo/parts/counterpart.py:

```
class CounterPart(Part):  
  
    counter = None # type: AttributeModel  
    delta = None  # type: AttributeModel  
  
    def zero(self):  
        """Zero the counter attribute"""  
        self.counter.set_value(0)  
  
    def increment(self):  
        """Add delta to the counter attribute"""  
        self.counter.set_value(self.counter.value +  
                                self.delta.value)
```



- Define the attributes as type *AttributeModel*
- Use the *AttributeModel*'s *set_value* method to update the counter's value
- Use its *value* attribute to get its value



Counter Part Definition

demo/parts/counterpart.py (continued):

```
def setup(self, registrar):  
    # type: (PartRegistrar) -> None  
    super(CounterPart, self).setup(registrar)  
    self.counter = NumberMeta(  
        "float64", "The current counter value",  
        tags=[config_tag(),  
              Widget.TEXTINPUT.tag()])  
    ).create_attribute_model()  
  
    self.delta = NumberMeta(  
        "float64", "Amount to increment by",  
        tags=[config_tag(),  
              Widget.TEXTINPUT.tag()])  
    ).create_attribute_model(  
        initial_value=1)
```

- Call the parent class *setup* method
- Create a *NumberMeta* object:
 - Type, description
 - Tags indicating usage info
- Create the *AttributeModel* from the meta object, supplying a default initial value if desired



Counter Part Definition

demo/parts/counterpart.py (continued):


```
registrar.add_attribute_model(  
    "counter",  
    self.counter,  
    self.counter.set_value)  
  
registrar.add_attribute_model(  
    "delta",  
    self.delta,  
    self.delta.set_value)  
  
registrar.add_method_model(self.zero)  
registrar.add_method_model(self.increment)
```

- Register the attribute models:
 - Name,
 - Model,
 - Method name
 - If not supplied, the attribute is read-only

- Register the method models as before



Running the Example





Practical Exercises



1. Modify the *Hello* example to throw an error saying “No name supplied!” if the name is an empty string.
2. Try setting the counter attribute to a non numeric value. What happens?
3. Modify the *Counter* example to add a new *decrement* method which decrements the counter by the delta value.