



Hardware Triggered Scanning: Introduction

Philip Taylor, Emma Arandjelovic
Observatory Sciences Limited




Course Aims

1. Get to know the functionality provided by the hardware triggered scanning stack
2. Understand the architecture, and acquire basic knowledge of all the components
3. Get hands-on experience of setting up scans
4. Learn how to debug common problems!



Course Content

1. Introduction and demo
2. Overview of the hardware
3. Low level control and data acq. concepts
4. Scan configuration
5. Experiment control and data processing
6. Practical exercises 




Grid Scanning Overview

- Move sample in X and Y to scan it through the beam
- Rewind at the end of each row or reverse direction ('alternate/snake' scan)

Data acquisition can be done using:

1. Software step scan
2. Hardware step scan
3. Continuous scan


Continuous scans are much faster! But also complex...





Framework Motivation

- Growing need for fast, continuous scanning
- Complex trajectories – not just grids!




- Need for on-the-fly visualization of data
- Desire for a generic solution



Hardware Components

Motion trajectory
control







Hardware Components

Motion trajectory
control

Data capture




Hardware Components

Motion trajectory
control

Flexible triggering, and
fast position capture

Data capture




Hardware Components

Motion trajectory
control

Flexible triggering, and
fast position capture

Data capture





Software Components



Data Analysis WorkbeNch

- Analysis and visualization




Generic Data Acquisition

- Experiment setup and supervision



Malcolm

- Scan configuration




Experimental Physics & Industrial Control System

- Low level control of hardware



System Architecture


Control Flow






System Architecture

Data Flow



Test Rig






Test Rig

- Beamline in a box!
 - Light illuminated 'sample'
 - Two stepper motors (θ and X)
 - DT Turbo Clipper controller
 - Allied Vision Mako camera
 - PandA electronics
 - PC with full software stack for testing and training





Test Rig






Demo Application

- Grid scan in θ and X
- 'Snake' style



Camera View




- Camera configured to capture a thin slice – each image is 1936 x 20 array
- Images used to reconstruct the sample

Stage 1: Setup scan




Open GDA
mapping
perspective



rotation axis

translation axis



Mapping Experiment Setup

Configure Beamline

Other Scan Axes

☐ stagez

Detectors

☐ Software-triggered camera 0.1

☒ Malcolm real motors 0.008

Region shape: Rectangle Scan path: Grid

stagea Start 0 mm stagea Points 60

stagea Stop 360 mm stagea Points 20

stagex Start -20 mm

stagex Stop 20 mm

Alternating ☒

Continuous ☒

Orientation Horizontal

Open DET panel from the Test Rig panel and select MJPG tab



Stage 2: Start scan

- Hit 'Queue scan'
- Motors move to initial positions
- Scan starts
- Images are acquired
- Watch the data acq. in GDA


Jython Console Queue Live Controls

Name	Status	Complete	Date Submitted	Mess
Sample 2 - Grid Scan	RUNNING	20%	11-Dec-2018 11:42:41	Point
Sample 2 - Grid Scan	COMPLETE	100%	10-Dec-2018 16:28:57	Scan
Sample 2 - Grid Scan	COMPLETE	100%	06-Dec-2018 15:22:22	Scan





Stage 3: Check the data






To improve the image:

- Higher resolution (AKA more images)
 - Increase stagea points → 180
 - Increase stagex points → 40
- Process the image
- If you see “jaggy edges” on the sinogram, try turning off alternating

Stage 4: Check the data





Stage 5: Add processing

Mapping Experiment Setup

Configure Beamline

Other Scan Axes

☐ stagez

Detectors

☐ Software-triggered camera 0.1

☒ Malcolm real motors 0.008

Region shape: Rectangle Scan path: Grid

stagea Start 0 mm stagea Points 60 - +

stagea Stop 360 mm stageex Points 20 - +

stagex Start -20 mm Alternating ☒

stagex Stop 20 mm Continuous ☒

Orientation Horizontal

Sample Name Edit metadata...

Processing

Processing Template and Detector Selection

Select the processing template file to use and the detector to apply it to.

Detector: Software-triggered camera

☐ Create a new processing file from a template:

Processing Template File: mean_integrate.nxs

☒ Use an existing processing file:

Processing File: /dls_sw/htss/software/gda

☐ Specify application and config. file:

App Name:

Config File:


Processing

App	Name
<input checked="" type="checkbox"/> dawn	mean_integrate.nxs

/dls_sw/htss/software/gda_var/processingTemplates/mean_integrate.nxs



Stage 6: Check the data






Stage 7: visualhulls

Processing Template and Detector Selection
Select the processing template file to use and the detector to apply it to.

Detector:

☐ Create a new processing file from a template:

Processing Template File: 

☐ Use an existing processing file:

Processing File:

☒ Specify application and config. file:

App Name:

Config File:


Processing		<input data-bbox="1593 535 1767 571" type="button" value="Add Processing..."/>
App	Name	
<input checked="" type="checkbox"/>	dawn	mean_integrate.nxs
<input checked="" type="checkbox"/>	visualhulls	vh.json

App Name must be **exactly** “visualhulls”
Config File: \$HOME/vh.json




Stage 7: visualhulls

1. Right click on the visualhulls .nxs in *Mapped Data*
2. Transfer > DataVis
3. Check the visualhulls.nxs file just created in *Data Files*
4. Select /processed/result and Plot Type “Volume” in *Datasets*





Stage 8: Check the data




Stage 9: Improve visualhulls

```

1 {
2     "step": 10,
3     "start": 400,
4     "stop": 1390,
5     "threshold": 230,
6     "cor": 0.8
7 }

```

Tab Width: 8 Ln 1, Col 1 INS



$\text{height}/2 = \text{cor}$

